



# Technical Alert

## Medical Device Recall: BIOFIRE FILMARRAY Pneumonia Panel

**TO:** Medical Staff and Clients

**FROM:** Dr. Amy Woron, PhD, MPH  
V.P. - Technical Director (Microbiology & Molecular Labs)

Jason Pon MB(ASCP), CHS(ABHI)  
Manager, DLS Molecular Labs

Dr. Wesley Kim, MD  
Medical Director DLS  
and QMC West

Dr. Ana Ortega-Lopez, MD  
Medical Director QMC Punchbowl,  
North Hawaii and Molokai

**DATE:** December 21, 2023

**SUBJECT: Medical Device Recall: False Positive Coronavirus on Biofire Filmarray Pneumonia Panel**

This alert applies to DLS ordering code #7212, Pneumonia Panel by Filmarray, which is only available to Emergency Departments and Inpatients.

BioMerieux, the manufacturer of BioFire FilmArray multiplex panels has issued an urgent medical device recall of the Pneumonia Panel because of an increased risk of false positive **seasonal** coronavirus results. This panel does not test for SARS-CoV-2 (COVID-19).

**IMPORTANT:**

**URGENT: MEDICAL DEVICE RECALL**

BIOFIRE® FILMARRAY® Pneumonia Panel – Ref. Numbers RFIT-ASY-0144 (30-pack) & RFIT-ASY-0145 (6-pack)

BIOFIRE® FILMARRAY® Pneumonia (PNplus) Panel plus – Ref. Numbers RFIT-ASY-0143 (30-pack) & RFIT-ASY-0142 (6 pack)

Increased risk of false positive seasonal Coronavirus results using the BIOFIRE® PN & PNplus Panel

December 2023

In the recall notice, bioMerieux reported that investigations indicated that the false positives appear to be as a result of non-specific amplification or cross-reactivity with high concentrations of human genomic DNA (hgDNA) that could be present in clinical specimens.

Although the notice indicates a recall, it does not appear that the manufacturer is pulling any product off the shelves. Consequently, **we are still able to offer the Pneumonia Panel** with some reporting restrictions. The manufacturer has advised customers that positive results for seasonal coronavirus should NOT BE REPORTED unless they can be confirmed. There is no confirmatory testing available, so positive results will not be reported. Negative results for seasonal coronavirus and all other pathogen targets will be reported as usual.

We do not have a timeframe in which the manufacturer expects to resolve this problem. We plan to update this alert as more information becomes available.

Please refer any questions to Jason Pon, Manager - DLS Molecular Laboratories at 808-441-5469 or DLS Client Services at 589-5101.

**Diagnostic Laboratory Services Inc.**

99-859 Iwaiwa Street – Aiea, HI 96701  
(808) 589-5100 – www.dlslab.com